

SGWO-45

High Concentration Rack Mounted Ozone Generator

KEY FEATURES

- Rack mounted
- 45G/HR ozone Output at 16LPM O₂
- Double 316 Stainless Steel and Quartz Electrodes Made in USA
- Tuned for Oxygen Only
- Built-in Flow Meter
- Power: 450W
- Weight: 9.7 KG
- CSA, ETL, UL
- 3 Year Warranty



PRODUCT DESCRIPTION

The ProMedUSA model SGWO-45 rack mounted ozone generator can produce up to 46 gram/hr of highly concentrated ozone. This unit is tuned for and designed to be used with an Oxygen (O₂) feed only.

The SGWO-45 utilizes our own 316 Stainless Steel and double Quartz technology electrodes (Made in USA) for ultra-pure and very highly concentrated ozone production

The SGWO-45 can be rack mounted, and requires the use of a Controller (SGWOC-R) and a balanced barometer (BBM).



Rack mounted 4 SGWO-45 Ozone Generators + 1 SGWOC-R Controller

SGWO-45 Ozone Production Graph

