

PROVIDER FEATURES

- High output pressure for medical applications
- Optional oxygen monitor
- Easy to service and maintain
- Time-proven NewLife® ruggedness and reliability
- Dual flow option available
- Dual function hour meter

PATIENT FEATURES

- Direct humidifier bottle attachment
- Resettable circuit breaker
- Simple user-interface
- Long tubing runs up to 200 ft (61 m)
- Low power draw 8 LPM model
- Backlined flow meter for easy viewing



ORDERING INFORMATION

Intensity Includes: Stationary Concentrator, DISS Fitting Adaptor, User Manual

INTENSITY 8

- AS094-100** 120V / 60Hz, USA
- AS094-101** 120V / 60Hz, USA, O₂ Monitor
- AS094-1** Dual Flow, 120V / 60Hz, USA
- AS094-4** Dual Flow, 120V / 60Hz, USA, O₂ Monitor

- AS094-103** 230V / 50 Hz, EU Cord, O₂ Monitor
- AS094-210** 230V / 50 Hz, UK Cord, O₂ Monitor
- AS094-10** Dual Flow, 230V / 50Hz, EU Cord, O₂ Monitor
- AS094-201** Dual Flow, 230V / 50Hz, UK Cord, O₂ Monitor

INTENSITY 10

- AS099-100** 120V / 60Hz, USA
- AS099-101** 120V / 60Hz, USA, O₂ Monitor
- AS099-1** Dual Flow, 120V / 60Hz, USA
- AS099-4** Dual Flow, 120V / 60Hz, USA, O₂ Monitor

- AS099-103** 230V / 50 Hz, EU Cord, O₂ Monitor
- AS099-210** 230V / 50 Hz, UK Cord, O₂ Monitor
- AS099-10** Dual Flow, 230V / 50Hz, EU Cord, O₂ Monitor
- AS099-201** Dual Flow, 230V / 50Hz, UK Cord, O₂ Monitor

*other models and country power cords available – [contact local customer service for details](#)

INTENSITY 8 SPECIFICATIONS

Flow Rates	2–8 LPM
O₂ Concentration	2–7 LPM: 92%+3.5%/-3%, 8 LPM: 90%+5.5%/-3%
Dimensions	27.5 x 16.5 x 14.5 in (69.9 x 41.9 x 36.8 cm)
Output Pressure	20 psig (138 kPA)
Weight	54 lbs (24.5 kg)
Power Consumption	410 W
Sound Level	52 dB(A)
Operating Environment	41°F to 104°F (5°C to 40°C), up to 95% humidity
Storage Environment	-4°F to 140°F (-20°C to 60°C), up to 95% relative humidity
Warranty	3 Years
Maintenance Schedule	Felt Filter – 1 Year Replacement, Intake Filter – Clean Weekly
Max Tubing	200 ft (61 m)

INTENSITY 10 SPECIFICATIONS

Flow Rates	2–10 LPM
O₂ Concentration	2–9 LPM: 92%+3.5%/-3%, 10 LPM: 90%+5.5%/-3%
Dimensions	27.5 x 16.5 x 14.5 in (69.9 x 41.9 x 36.8 cm)
Output Pressure	20 psig (138 kPA)
Weight	58 lbs (26.3 kg)
Power Consumption	590 W
Sound Level	55 dB(A)
Operating Environment	41°F to 104°F (5°C to 40°C), up to 95% humidity
Storage Environment	-4°F to 140°F (-20°C to 60°C), up to 95% relative humidity
Warranty	3 Years
Maintenance Schedule	Felt Filter – 1 Year Replacement, Intake Filter – Clean Weekly
Max Tubing	200 ft (61 m)